

A Publication of Reliable Methods for the Preparation of Organic Compounds

Working with Hazardous Chemicals

The procedures in Organic Syntheses are intended for use only by persons with proper training in experimental organic chemistry. All hazardous materials should be handled using the standard procedures for work with chemicals described in references such as "Prudent Practices in the Laboratory" (The National Academies Press, Washington, D.C., 2011; the full accessed of charge text can be free at http://www.nap.edu/catalog.php?record_id=12654). All chemical waste should be disposed of in accordance with local regulations. For general guidelines for the management of chemical waste, see Chapter 8 of Prudent Practices.

In some articles in *Organic Syntheses*, chemical-specific hazards are highlighted in red "Caution Notes" within a procedure. It is important to recognize that the absence of a caution note does not imply that no significant hazards are associated with the chemicals involved in that procedure. Prior to performing a reaction, a thorough risk assessment should be carried out that includes a review of the potential hazards associated with each chemical and experimental operation on the scale that is planned for the procedure. Guidelines for carrying out a risk assessment and for analyzing the hazards associated with chemicals can be found in Chapter 4 of Prudent Practices.

The procedures described in *Organic Syntheses* are provided as published and are conducted at one's own risk. *Organic Syntheses, Inc.,* its Editors, and its Board of Directors do not warrant or guarantee the safety of individuals using these procedures and hereby disclaim any liability for any injuries or damages claimed to have resulted from or related in any way to the procedures herein.

These paragraphs were added in September 2014. The statements above do not supersede any specific hazard caution notes and safety instructions included in the procedure.

Organic Syntheses, Coll. Vol. 8, p.347 (1993); Vol. 66, p.8 (1988).

A GENERAL [3+2] ANNULATION: cis-4-exo-ISOPROPENYL-1,9-DIMETHYL-8-(TRIMETHYLSILYL)BICYCLO[4.3.0]NON-8-EN-2-ONE

[4*H*-Inden-4-one, 1,3a,5,6,7,7a-hexahydro-3,3a-dimethyl-6-(1-methylethenyl)-2-(trimethylsilyl)-, (3aα,6α,7aα)-]



Submitted by Rick L. Danheiser, David M. Fink, and Yeun-Min Tsai¹. Checked by Marianne Marsi and Bruce E. Smart.

1. Procedure

A 500-mL, three-necked, round-bottomed flask is equipped with a 25-mL pressure-equalizing dropping funnel, a mechanical stirrer, and a Claisen adapter fitted with a nitrogen inlet adapter and a low-temperature thermometer (Note 1). The flask is charged with 115 g (0.077 mol) of (*R*)-(–)-carvone (Note 2), 10.8 g (0.079 mol) of 1-methyl-1-(trimethylsilyl)allene (Note 3), and 180 mL of dry dichloromethane (Note 4), and then cooled below -75° C with a dry ice–acetone bath while a solution of 17.4 g (0.092 mol) of titanium tetrachloride (Note 5) in 10 mL of dichloromethane is added dropwise over 1 hr. After 30 min, the cold bath is removed, and the reaction mixture, which appears as a red suspension, is allowed to warm to 0°C over approximately 30 min. The resulting dark red solution is poured slowly into a 2-L Erlenmeyer flask containing a magnetically stirred mixture of 400 mL of diethyl ether and 400 mL of water (Note 6). The aqueous phase is separated and extracted with two 200-mL portions of ether. The combined organic phases are washed with 250 mL of water and 250 mL of saturated sodium chloride solution, dried over anhydrous magnesium sulfate, and concentrated at reduced pressure using a rotary evaporator. The residual yellow liquid is distilled through a 15-cm Vigreux column at reduced pressure to afford 17.5 g (82%) of the bicyclonoenone **1** as a very-pale-yellow liquid, bp 98–101°C (0.03 mm), $[\alpha]_D^{20} -157.8 \pm 0.8 (1.57, CH_2Cl_2)$ (Note 7) and (Note 8).

2. Notes

1. The apparatus is flame-dried under vacuum and then maintained under an atmosphere of nitrogen during the course of the reaction.

2. (R)-(-)-Carvone was purchased from Aldrich Chemical Company, Inc. and distilled before use.

3. 1-Methyl-1-(trimethylsilyl)allene (90% purity, contaminated with 10% 1-trimethylsilyl-1-butyne) was prepared by the method of Danheiser, R. L.; Tasai, Y. M.; Fink, D. M. *Org. Synth., Coll. Vol. VIII*, **1993**, 471.

4. Dichloromethane was distilled from calcium hydride immediately before use.

5. Titanium tetrachloride (99.9%) was obtained from the Aldrich Chemical Company, Inc. and distilled before use. Lower yields (70–77%) resulted if the titanium tetrachloride was not distilled. Unreacted carvone is recovered if a small excess of titanium tetrachloride is not used.

6. The two-phase mixture is vigorously stirred using a 7-cm Teflon-coated magnetic stirring bar.

7. The submitters report obtaining 18.8 g (88%) of product, bp 93–96°C (0.03 mm). The purity of the product was determined to be >99% by gas chromatographic analysis (10% OV-101 on 100-120 mesh Chromosorb W, 6 ft × 1/8 in., program: 200°C for 2 min and then 200–300°C at 32°C/min).

8. The product exhibits the following spectral properties: IR (neat) cm⁻¹: 3080, 2950, 2920, 1700, 1640, 1610, 1440, 1375, 1315, 1245, 1200, 830, 755, 680; ¹H NMR (250 MHz, CDCl₃) δ: 0.08 (s, 9 H), 1.10

(s, 3 H), 1.65 (t, 3 H, J = 2.2), 1.68 (m, 3 H), 1.65–1.72 (m, 2 H), 2.10 (d of m, 1 H, J = 12.4), 2.15–2.25 (m, 1 H), 2.27–2.31 (m, 2 H), 2.45–2.57 (m, 2 H), 4.63 (m, 1 H), 4.72 (m, 1 H); ¹³C NMR (62.8 MHz, CDCl₃) δ : -0.7, 14.4, 21.1, 21.3, 32.4, 39.8, 42.0, 43.9, 46.1, 64.8, 110.4, 136.1, 147.7, 151.5, 215.6.

3. Discussion

The procedure described here serves to illustrate a general [3 + 2] annulation method^{2 3} for the synthesis of cyclopentane derivatives. A unique feature of this one-step annulation is its capacity to generate regiospecifically five-membered rings substituted at each position, and functionally equipped for further synthetic elaboration. As formulated in the following equation, the reaction proceeds with remarkably high stereoselectivity via the effective suprafacial addition of the three-carbon allene component to an electron-deficient olefin ("allenophile").



Some representative examples of the [3 + 2] annulation are listed in Table I. Both cyclic and acyclic allenophiles participate in the reaction. α -Alkylidene ketones undergo annulation to provide access to spiro-fused systems, and acetylenic allenophiles react to form cyclopentadiene derivatives. The reactions of (*E*)- and (*Z*)-3-methyl-3-penten-2-one illustrate the stereochemical course of the annulation, which proceeds with a strong preference for the suprafacial addition of the allene to the two-carbon allenophile. The high stereoselectivity displayed by the reaction permits the stereocontrolled synthesis of a variety of mono- and polycyclic systems.

L			
Allenophile	Allene	Annulation Product	Yield (%)
O H	l ₂ C=C=C ^{SIMe₃}	SiMe3	71-75
		SiMe3	86
MeO		H ₃ CO-C SiMe ₃	49
		SiMe3	91
,)			90
		SiMe3	71

TABLE I[3 + 2] ANNULATIONS EMPLOYING ALLENYLSILANES



This preparation is referenced from:

• Org. Syn. Coll. Vol. 8, 471

References and Notes

- 1. Department of Chemistry, Massachusetts Institute of Technology, Cambridge, MA 02139.
- 2. Danheiser, R. L.; Carini, D. J.; Basak, A. J. Am. Chem. Soc. 1981, 103, 1604;
- 3. Danheiser, R. L.; Carini, D. J.; Fink, D. M.; Basak, A. Tetrahedron 1983, 39, 935.

Appendix Chemical Abstracts Nomenclature (Collective Index Number); (Registry Number)

(R)-(-)-carvone

(E)- and (Z)-3-methyl-3-penten-2-one

ether,

diethyl ether (60-29-7)

sodium chloride (7647-14-5)

nitrogen (7727-37-9)

dichloromethane (75-09-2)

magnesium sulfate (7487-88-9)

Allene (463-49-0)

titanium tetrachloride (7550-45-0)

calcium hydride (7789-78-8)

carvone

1-Methyl-1-(trimethylsilyl)allene (74542-82-8)

1-trimethylsilyl-1-butyne

cis-4-exo-Isopropenyl-1,9-dimethyl-8-(trimethylsilyl)bicyclo[4.3.0]non-8-en-2-one (77494-23-6)

4H-Inden-4-one, 1,3a,5,6,7,7a-hexahydro-3,3a-dimethyl-6-(1-methylethenyl)-2-(trimethylsilyl)-, $(3a\alpha,6\alpha,7a\alpha)$ -

Copyright © 1921-2005, Organic Syntheses, Inc. All Rights Reserved