



A Publication
of Reliable Methods
for the Preparation
of Organic Compounds

Working with Hazardous Chemicals

The procedures in *Organic Syntheses* are intended for use only by persons with proper training in experimental organic chemistry. All hazardous materials should be handled using the standard procedures for work with chemicals described in references such as "Prudent Practices in the Laboratory" (The National Academies Press, Washington, D.C., 2011; the full text can be accessed free of charge at http://www.nap.edu/catalog.php?record_id=12654). All chemical waste should be disposed of in accordance with local regulations. For general guidelines for the management of chemical waste, see Chapter 8 of Prudent Practices.


In some articles in *Organic Syntheses*, chemical-specific hazards are highlighted in red "Caution Notes" within a procedure. It is important to recognize that the absence of a caution note does not imply that no significant hazards are associated with the chemicals involved in that procedure. Prior to performing a reaction, a thorough risk assessment should be carried out that includes a review of the potential hazards associated with each chemical and experimental operation on the scale that is planned for the procedure. Guidelines for carrying out a risk assessment and for analyzing the hazards associated with chemicals can be found in Chapter 4 of Prudent Practices.

The procedures described in *Organic Syntheses* are provided as published and are conducted at one's own risk. *Organic Syntheses, Inc.*, its Editors, and its Board of Directors do not warrant or guarantee the safety of individuals using these procedures and hereby disclaim any liability for any injuries or damages claimed to have resulted from or related in any way to the procedures herein.

September 2014: The paragraphs above replace the section "Handling and Disposal of Hazardous Chemicals" in the originally published version of this article. The statements above do not supersede any specific hazard caution notes and safety instructions included in the procedure.

Organic Syntheses, Coll. Vol. 9, p.472 (1998); Vol. 72, p.180 (1995).

PREPARATION AND REACTIONS OF ALKENYLCHROMIUM REAGENTS: 2-HEXYL-5-PHENYL-1-PENTEN-3-OL



Submitted by Kazuhiko Takai¹, Koichi Sakogawa, Yasutaka Kataoka, Koichiro Oshima, and Kiiro Utimoto².

Checked by Tadahiro Takemoto and Larry E. Overman.

1. Procedure

A. *1-Hexylethenyl triflate* (Note 1). A dry, 300-mL, three-necked, round-bottomed flask is equipped with a magnetic stirring bar, nitrogen inlet, rubber septum, and a 100-mL, graduated, pressure-equalizing addition funnel that is sealed with a rubber septum. The flask is charged with *1-octyne* (11 g, 0.10 mol, (Note 2)) and dry *hexane* (100 mL, (Note 3)). The contents of the flask are cooled to -30°C (dry ice-methanol). The addition funnel is charged with neat *trifluoromethanesulfonic acid* (TfOH, 4.9 mL, 55 mmol, (Note 4)) which is then added dropwise at -30°C to the stirring solution over a period of 30 min. The addition funnel is rinsed with dry *hexane* (10 mL). The dark brown solution is warmed to 0°C and stirred at 0°C for 10 min. Saturated *sodium bicarbonate* (NaHCO_3) solution (20 mL) is introduced at 0°C to the reaction mixture. The organic layer is separated, washed with additional NaHCO_3 solution (2×50 mL), and then dried over *potassium carbonate*. The solvent is removed with a rotary evaporator. The crude product is distilled using a short-path still to give 9.6–10.4 g (67–73%) of *1-hexylethenyl triflate* (Note 5), bp $97\text{--}98^{\circ}\text{C}$ at 19 mm.

B. *2-Hexyl-5-phenyl-1-penten-3-ol*. A dry, 500-mL, four-necked, round-bottomed flask is equipped with a mechanical stirring bar, nitrogen inlet, rubber septum, and a 100-mL, graduated, pressure-equalizing addition funnel that is sealed with a rubber septum. In the flask are placed anhydrous *chromium(II) chloride* (CrCl_2) (10 g, 80 mmol, (Note 6)) and anhydrous *nickel(II) chloride* (NiCl_2) (52 mg, 0.40 mmol, (Note 7)) under an *argon* atmosphere. The flask is cooled to 0°C and dry, oxygen-free *N,N*-dimethylformamide (DMF, 250 mL, (Note 8)) is added to the flask with stirring. The salts are dissolved in a slightly exothermic process. The mixture is stirred at 0°C for 10 min. To the $\text{CrCl}_2\text{-NiCl}_2$ reagent at 25°C is added a solution of *3-phenylpropanal* (2.7 g, 20 mmol, (Note 9)) in DMF (20 mL) by syringe. A solution of *1-hexylethenyl triflate* (10 g, 40 mmol, (Note 5)) in DMF (60 mL) is added at 25°C through the addition funnel over a period of 5 min. The entire mixture is stirred at 25°C for 30 min. The reaction mixture is diluted with *ether* (200 mL), poured into ice-cooled water (400 mL), and extracted with *ether* (3×200 mL) repeatedly. The combined extracts are washed with aqueous *sodium chloride* solution (150 mL), dried over anhydrous *sodium sulfate*, and concentrated with a rotary evaporator (25°C , water bath). The crude product is distilled using a short-path still to give 4.0–4.6 g (82–94%) of *2-hexyl-5-phenyl-1-penten-3-ol*, bp $109\text{--}111^{\circ}\text{C}$ at 0.11 mm (Note 10) and (Note 11).

2. Notes

1. This procedure was reported by Stang and Summerville.³
2. *1-Octyne* was distilled, bp $125\text{--}126^{\circ}\text{C}$.
3. *Hexane* was freshly distilled from *sodium* or *calcium hydride*.
4. *Trifluoromethanesulfonic acid* was distilled, bp $72\text{--}73^{\circ}\text{C}$ at 20 mm.

5. The distilled triflate contained 3–9% of isomers, that were confirmed by capillary GLPC (Silicone OV-17, 50 m, 105°C, **1-hexylethenyl triflate**: $t_R = 11.1$ min; isomers: $t_R = 11.4$ and 12.4 min) and ^1H NMR analysis. Because the isomers did not interfere with the second reaction, the triflate was employed without further purification. Spectral data of the distilled triflate was as follows: IR (neat) cm^{-1} : 2954, 2930, 2858, 1671, 1419, 1250, 1213, 1141, 943, 899, 703, 610; ^1H NMR (CDCl_3) δ : 0.89 (t, 3 H, $J = 6.9$), 1.2–1.4 (m, 6 H), 1.5–1.6 (m, 2 H), 2.34 (t, 2 H, $J = 7.5$), 4.93 (dt, 1 H, $J = 3.5, 1.0$), 5.08 (d, 1 H, $J = 3.5$), ^1H NMR peaks of the impurities appeared at δ 2.05–2.20 and 5.18–5.55; ^{13}C NMR (CDCl_3) δ : 13.9, 22.4, 25.9, 28.3, 31.3, 33.8, 103.9, 118.5 (q, $J = 320, \text{CF}_3$), 157.1.
6. **Chromium(II) chloride** (95% purity) was purchased from Aldrich Chemical Company, Inc., and used without further purification. The salt is easily oxidized and should be handled under an inert atmosphere.
7. Anhydrous **nickel(II) chloride** was purchased from Nacalai Tesque Co. and used without further purification. The salt is hygroscopic and should be handled under an inert atmosphere.
8. **N,N-Dimethylformamide** was refluxed in the presence of **calcium sulfate** under reduced pressure and distilled with **nitrogen** bubbling from a capillary, bp 76°C at 39 mm.
9. **3-Phenylpropanal** was freshly distilled, bp 97–98°C at 12 mm.
10. The pot residue can be bulb-to-bulb distilled to give an additional 0.5–1 g of product to bring the combined yields to 94–95%. **2-Hexyl-5-phenyl-1-penten-3-ol** has the following properties: $R_f = 0.29$ (**ethyl acetate/hexane** = 1/10); IR (neat) cm^{-1} : 3340, 3024, 2924, 2854, 1647, 1600, 1497, 1456, 1017, 900, 741, 697; ^1H NMR (CDCl_3) δ : 0.88 (t, 3 H, $J = 6.8$), 1.2–1.4 (m, 8 H), 1.7–2.2 (m, 5 H), 2.62 (ddd, 1 H, $J = 6.6, 9.6, 13.9$), 2.74 (ddd, 1 H, $J = 6.1, 9.6, 13.9$), 4.10 (dd, 1 H, $J = 5.1, 7.5$), 4.87 (d, 1 H, $J = 1.5$), 5.04 (s 1 H), 7.1–7.4 (m, 5 H). ^{13}C NMR (CDCl_3) δ : 14.1, 22.6, 27.9, 29.2, 31.4, 31.7, 31.9, 37.0, 74.7, 109.3, 125.7, 128.3, 128.4, 142.0, 152.0.
11. The following ratio of reactants, aldehyde/alkenyl triflate/ $\text{CrCl}_2/\text{NiCl}_2 = 1/2/4/0.02$ gave the best results. When the ratio of reagents was reduced to aldehyde/alkenyl triflate/ $\text{CrCl}_2/\text{NiCl}_2 = 1/1/2/0.01$, the reaction proceeded slowly. When the reaction was carried out at 25°C for 2 hr, 65% of the **2-hexyl-5-phenyl-1-penten-3-ol** was isolated and 14% (GLPC analysis) of the **3-phenylpropanal** remained.

Waste Disposal Information

Compounds of chromium and nickel are toxic. The aqueous layers from Step B and any other waste materials should be disposed of properly (see "Prudent Practices in the Laboratory"; National Academy Press; Washington, DC, 1995).

3. Discussion




















In the synthesis of a complex molecule, it is sometimes necessary to prepare an organometallic reagent under mild conditions. Because of the strong basicity and nucleophilicity of alkenyllithium and alkenylmagnesium compounds, only a few electrophilic functional groups are stable under the reaction conditions.⁴ In contrast, the alkenylchromium reagents described here react with aldehydes to give adducts in good to excellent yields in the presence of ketones, esters, amides, acetals, ethers, silyl ethers, and nitriles.^{5,6,7} The method is especially effective for highly oxygenated molecules.⁶ Intramolecular cyclization of iodo aldehydes leading to 13-membered lactones has also been reported.⁸


Addition of a catalytic amount of NiCl_2 to CrCl_2 is essential for the formation of alkenylchromium reagents.^{5,6,7} However, a substantial amount of 1,3-diene, the coupling product of the alkenyl iodide, is produced if a higher content of NiCl_2 is employed.^{9,10}

Reduction of alkenyl iodides and bromides to alkenylchromium reagents with CrCl_2 proceeds smoothly under the same conditions.^{5,6,7} Several examples of the Grignard-type addition of alkenyl halides and triflates to aldehydes with the combination of CrCl_2 and NiCl_2 are shown in the Table. Iodoalkenes are more reactive than bromoalkenes (compare run 1 and 2). In the case of an α,β -unsaturated aldehyde, the 1,2-addition product is the main product (run 4). The alkenylchromium reagent adds to an aldehyde group selectively (runs 5–7), as do allyl-¹¹ and alkynylchromium reagents.¹² Steric interaction of substituents at a position cis to halogen causes cis-trans isomerization in some cases.¹⁰ For example, while the reaction of (E)- and (Z)-2-bromostyrene and **benzaldehyde** proceeded stereospecifically (runs 8 and 9), both (E)- and (Z)-2-iodo-1-phenyl-1-propene reacted with

benzaldehyde to give (E)-1,3-diphenyl-2-methyl-2-propen-1-ol as the sole product (runs 10 and 11). The regiochemistry of double bonds is maintained during the coupling reaction. The CrCl₂-NiCl₂ system is also effective for the addition of iodobenzene to an aldehyde (run 12).⁷

TABLE
GRIGNARD-TYPE COUPLING BETWEEN ALKENYL HALIDES (OR TRIFLATES) AND ALDEHYDES, MEDIATED BY THE CrCl₂-NiCl₂ SYSTEM^a

Run	Alkenyl Halide (or Triflate)	Aldehyde	Time hr	Product	Yield %
				$\text{R}^1\text{C}(\text{X})=\text{C} + \text{R}^2\text{CHO} \xrightarrow[\text{DMF, 25}^\circ\text{C}]{\text{CrCl}_2, \text{cat. NiCl}_2} \text{R}^1\text{C}(\text{OH})\text{C}(\text{R}^2)=\text{C}$	
1		n-C ₈ H ₁₇ CHO	0.25		100
2		n-C ₈ H ₁₇ CHO	3		62
3		n-C ₈ H ₁₇ CHO	3		81
4		(E)-PrCH=CHO	4		64
5			1		87
6			2		78
7		PhCOMe	6	(recovery of PhCOMe, 87%)	
8		PhCHO	1		82
9		PhCHO	1		78
10		PhCHO	3		91 ^d
11		PhCHO	3		90 ^d
12	PhI	n-C ₈ H ₁₇ CHO	3		83



^aA mixture of an alkenyl halide (or alkenyl triflate, 2.0 mmol) and an aldehyde (1.0 mmol) was treated at 25°C with CrCl₂ (4.0 mmol) and NiCl₂ (0.02 mmol) in DMF.

^bProducts are isolated by column chromatography on silica gel. ^cThe reaction mixture was heated at 60°C. ^dThe absence of (Z)-isomer was confirmed by GLPC (Silicone SE-30) and NMR analysis.

References and Notes

1. Department of Applied Chemistry, Faculty of Engineering, Okayama University, Tsushima, Okayama 700, Japan.
2. Division of Material Chemistry, Faculty of Engineering, Kyoto University, Yoshida, Kyoto 606-01, Japan.
3. Stang, P. J.; Summerville, R. H. *J. Am. Chem. Soc.* **1969**, *91*, 4600; Summerville, R. H.; Schleyer, P. v. R. *J. Am. Chem. Soc.* **1974**, *96*, 1110.
4. Recently, functionalized alkenyllithiums could be prepared from the corresponding alkenyl iodides by transmetalation in THF:ether:pentane (4:1:1) at -100°C, Tucker, C. E.; Majid, T. N.; Knochel, P. *J. Am. Chem. Soc.* **1992**, *114*, 3983.
5. Takai, K.; Tagashira, M.; Kuroda, T.; Oshima, K.; Utimoto, K.; Nozaki, H. *J. Am. Chem. Soc.* **1986**, *108*, 6048.
6. Jin, H.; Uenishi, J.-i.; Christ, W. J.; Kishi, Y. *J. Am. Chem. Soc.* **1986**, *108*, 5644; Aicher, T. D.; Buszek, K. R.; Fang, F. G.; Forsyth, C. J.; Jung, S. H.; Kishi, Y.; Matelich, M. C.; Scola, P. M.; Spero, D. M.; Yoon, S. K. *J. Am. Chem. Soc.* **1992**, *114*, 3162.
7. Takai, K.; Kimura, K.; Kuroda, T.; Hiyama, T.; Nozaki, H. *Tetrahedron Lett.* **1983**, *24*, 5281. In 1983 the authors/submitters purchased anhydrous CrCl₂ from ROC/RIC Corp. (507-519 Main St., Belleville, NJ 07109) and used it without further purification. The effective lots contained ca. 0.5 mol% of Ni based on Cr.
8. Schreiber, S. L.; Meyers, H. V. *J. Am. Chem. Soc.* **1988**, *110*, 5198.
9. Kende, A. S.; Liebeskind, L. S.; Braitsch, D. M. *Tetrahedron Lett.* **1975**, 3375; Zembayashi, M.; Tamao, K.; Yoshida, J.-i.; Kumada, M. *Tetrahedron Lett.* **1977**, 4089.
10. The cis-trans isomerization of the double bond could occur during oxidative addition of triflate to nickel(0), Semmelhack, M. F.; Helquist, P. M.; Gorzynski, J. D. *J. Am. Chem. Soc.* **1972**, *94*, 9234.
11. Hiyama, T.; Okude, Y.; Kimura, K.; Nozaki, H. *Bull. Chem. Soc. Jpn.* **1982**, *55*, 561 and references cited therein.
12. Takai, K.; Kuroda, K.; Nakatsukasa, S.; Oshima, K.; Nozaki, H. *Tetrahedron Lett.* **1985**, *26*, 5585.

Appendix Chemical Abstracts Nomenclature (Collective Index Number); (Registry Number)

oxygen-free N,N-dimethylformamide


CrCl₂-NiCl₂

NiCl₂



(E)- and (Z)-2-bromostyrene

(E)- and (Z)-2-iodo-1-phenyl-1-propene



Phl

potassium carbonate (584-08-7)

ethyl acetate (141-78-6)

ether (60-29-7)

sodium bicarbonate,
 NaHCO_3 (144-55-8)

sodium chloride (7647-14-5)

sodium sulfate (7757-82-6)

nitrogen (7727-37-9)

benzaldehyde,
 PhCHO (100-52-7)

calcium sulfate (7778-18-9)

nickel(0) (7440-02-0)

sodium (13966-32-0)

Iodobenzene (591-50-4)

N,N-dimethylformamide,
DMF (68-12-2)

hexane (110-54-3)

chromium(II) chloride (10049-05-5)

argon (7440-37-1)

calcium hydride (7789-78-8)

nickel(II) chloride (7718-54-9)

1-OCTYNE (629-05-0)

triflate

trifluoromethanesulfonic acid (1493-13-6)

3-phenylpropanal (104-53-0)

2-Hexyl-5-phenyl-1-penten-3-ol (187821-45-0)

1-Hexylethenyl triflate (98747-02-5)

(E)-1,3-diphenyl-2-methyl-2-propen-1-ol