



A Publication
of Reliable Methods
for the Preparation
of Organic Compounds

Working with Hazardous Chemicals

The procedures in *Organic Syntheses* are intended for use only by persons with proper training in experimental organic chemistry. All hazardous materials should be handled using the standard procedures for work with chemicals described in references such as "Prudent Practices in the Laboratory" (The National Academies Press, Washington, D.C., 2011; the full text can be accessed free of charge at http://www.nap.edu/catalog.php?record_id=12654). All chemical waste should be disposed of in accordance with local regulations. For general guidelines for the management of chemical waste, see Chapter 8 of Prudent Practices.

In some articles in *Organic Syntheses*, chemical-specific hazards are highlighted in red "Caution Notes" within a procedure. It is important to recognize that the absence of a caution note does not imply that no significant hazards are associated with the chemicals involved in that procedure. Prior to performing a reaction, a thorough risk assessment should be carried out that includes a review of the potential hazards associated with each chemical and experimental operation on the scale that is planned for the procedure. Guidelines for carrying out a risk assessment and for analyzing the hazards associated with chemicals can be found in Chapter 4 of Prudent Practices.


The procedures described in *Organic Syntheses* are provided as published and are conducted at one's own risk. *Organic Syntheses, Inc.*, its Editors, and its Board of Directors do not warrant or guarantee the safety of individuals using these procedures and hereby disclaim any liability for any injuries or damages claimed to have resulted from or related in any way to the procedures herein.

September 2014: The paragraphs above replace the section "Handling and Disposal of Hazardous Chemicals" in the originally published version of this article. The statements above do not supersede any specific hazard caution notes and safety instructions included in the procedure.

Copyright © 2007 Organic Syntheses, Inc. All Rights Reserved

**GENERATION OF YNOLATE AND Z-SELECTIVE OLEFINATION
OF ACYLSILANES: (Z)-2-METHYL-3-TRIMETHYLSILYL-2-
BUTENOIC ACID**

((Z)-2-Butenoic acid, 2-methyl-3-trimethylsilyl-)



Submitted by Mitsuru Shindo,¹ Kenji Matsumoto,² and Kozo Shishido.²

Checked by Scott E. Denmark and Brian M. Eklov.

1. Procedure

Caution! tert-Butyllithium is extremely pyrophoric and must not be allowed to come into contact with the atmosphere. This reagent should only be handled by individuals trained in its proper and safe use. It is recommended that transfers be carried out by using a 20-mL or smaller glass syringe filled to no more than 2/3 capacity, or by cannula. For a discussion of procedures for handling air-sensitive reagents, see Aldrich Technical Bulletin AL-134. [Note added August 2009].

A. *Ethyl 2,2-dibromopropionate (1).*³ A 1-L, three-necked, round-bottomed flask is equipped with a magnetic stirring bar and fitted with a thermocouple temperature probe (Note 1) inserted through an adapter, a rubber septum, and a 100-mL, pressure-equalizing addition funnel fitted with a vacuum adapter which is capped with a rubber septum and the sidearm of which is connected to a vacuum/argon gas line. The flask is flame-dried under reduced pressure and then is maintained under an

atmosphere of argon during the course of the reaction. The flask is charged with diisopropylamine (23.0 mL, 164 mmol, 1.1 equiv) (Note 2) and anhydrous THF (350 mL) (Note 3), and then is cooled to $-71\text{ }^{\circ}\text{C}$ (internal) in a dry ice/isopropyl alcohol bath (Note 4) while a 1.62 M solution of *n*-butyllithium in hexane (102 mL, 165 mmol, 1.1 equiv) (Note 5) is added through the addition funnel over 10 min. After an additional 16 min, a solution of ethyl 2-bromopropanoate (20.0 mL, 27.8 g, 154 mmol) (Note 6) in anhydrous THF (20 mL) is added dropwise through the addition funnel over 18 min (the addition funnel is rinsed with 10 mL of anhydrous THF). After an additional 34 min, a solution of 1,2-dibromo-1,1,2,2-tetrachloroethane (53.7 g, 165 mmol, 1.1 equiv) (Note 7) in THF (60 mL) is added rapidly in one portion through the addition funnel with vigorous stirring (Note 8). After being stirred in the cold bath for 40 min (Note 9), the mixture is poured into 150 mL of saturated aqueous sodium bicarbonate solution in a 1-L separatory funnel (Note 10) and the mixture is extracted with hexanes (2 x 150 mL). The combined organic layers are washed with 150 mL of water and 150 mL of brine, then are dried over MgSO_4 (6 g), filtered under vacuum, and concentrated by rotary evaporation ($40\text{ }^{\circ}\text{C}$, 50–100 mmHg) to afford a dark brown oil. The oil is transferred to a 100-mL, round-bottomed flask equipped with a magnetic stirbar and fitted with a short-path distillation unit. Distillation of the product under reduced pressure using an oil bath affords 36.1 g (90% yield) of ethyl 2,2-dibromopropionate (**1**) as a colorless liquid ($71\text{--}80\text{ }^{\circ}\text{C}$ / 17 mmHg) (Notes 11, 12) after an approximate 5 g forerun.

*B. (Z)-2-Methyl-3-trimethylsilyl-2-butenic acid (2).*⁴ A 500-mL, three-necked, round-bottomed flask is equipped with a magnetic stirring bar and fitted with a thermocouple temperature probe adapter, a rubber septum, and a pressure-equalizing, 200-mL addition funnel fitted with a vacuum adapter capped with a rubber septum. The side arm of the vacuum adapter is connected to a vacuum/argon gas system. The flask is flame-dried under reduced pressure and then maintained under an atmosphere of argon during the course of the reaction. The thermocouple is inserted through the adapter and the flask is charged with ethyl 2,2-dibromopropionate (**1**, 17.6 g, 68.0 mmol, 1.2 equiv) and anhydrous THF (180 mL) (Note 3). The solution is

cooled with stirring to $-72\text{ }^{\circ}\text{C}$ (internal) in a dry ice/isopropyl alcohol bath, then a 1.72 M solution of *tert*-butyllithium in pentane (158 mL, 272 mmol, 4.9 equiv) (Note 13) is added dropwise over 67 min through the addition funnel such that the internal temperature is kept below $-64\text{ }^{\circ}\text{C}$ (Note 14). The addition funnel is rinsed with 5 mL of THF. The yellow solution is stirred at $-72\text{ }^{\circ}\text{C}$ for an additional 110 min (3 h total reaction time), and then the cooling bath is replaced by an ice bath. After being stirred at $2\text{ }^{\circ}\text{C}$ (internal) for 30 min, the ice bath is removed, and the mixture is allowed to warm to room temperature. A solution of acetyltrimethylsilane (6.50 g, 56 mmol, Note 15) in anhydrous THF (15 mL) is added through the addition funnel over 15 min (Note 16). The addition funnel is rinsed with 5 mL of anhydrous THF, and the mixture is stirred for 45 min at room temperature. The reaction mixture is poured onto 1 M aqueous NaOH (120 mL) in a 1-L separatory funnel. The organic phase is separated and the aqueous phase (pH 8) is washed with ethyl acetate (2 x 120 mL). The combined organic extracts are washed with 1 M aqueous sodium hydroxide solution (120 mL) and set aside. The combined aqueous extracts (pH 9) are acidified with 6 M HCl (100 mL) to pH 1. The resulting cloudy suspension is extracted with ethyl acetate (3 x 120 mL). These organic extracts are combined, washed with brine (120 mL), dried over MgSO_4 (4 g), filtered through a coarse glass frit and then are concentrated on a rotary evaporator ($23\text{ }^{\circ}\text{C}$, 10 mmHg) to yield 7.64 g of a viscous yellow-orange oil that solidifies upon standing in the freezer ($-20\text{ }^{\circ}\text{C}$) (Note 17). The crude product is transferred to a sublimation chamber (chamber: inside diameter, 57 mm; height, 90 mm; condenser: outside diameter, 27 mm; height 70 mm). The bottom of the chamber is cooled in an ice bath and the apparatus is attached to a vacuum line. The bath is removed once the majority of the volatile materials are evaporated as judged by stabilization of the pressure. When the pressure inside the chamber drops below 0.1 mmHg, the condenser is charged with a dry ice/isopropyl alcohol mixture and the chamber is placed in an oil bath, pre-heated to $40\text{ }^{\circ}\text{C}$. The sublimate was removed from the cold finger three times to afford 6.12 – 6.60 g (64–69%) of **2** as a white crystalline solid (mp $48.5\text{--}50.5\text{ }^{\circ}\text{C}$, Note 19).

2. Notes

1. A PFA-coated thermocouple probe, Type K (Omega Engineering, Inc.) was inserted through the septum.

2. Diisopropylamine was purchased from Sigma-Aldrich Company (>98.0%) and was distilled from calcium hydride.

3. Anhydrous THF (>99.5%, stabilizer free) was purchased from Fischer Scientific and was used as received (Karl Fischer analysis gave < 120 µg/mL).

4. For a small scale (smaller than 500-mL flask), the submitters used a low temperature bath with magnetic stirrer (Tokyo Rikakikai Co., Ltd., PSL-1800).

5. *n*-Butyllithium solution was purchased from Sigma-Aldrich and was titrated by both No-D NMR spectroscopy^{5a} and Gilman double titration.^{5b}

6. Ethyl 2-bromopropanoate (>98%) was purchased from Tokyo Kasei Kogyo Co. Ltd. (TCI) and was distilled.

7. 1,2-Dibromo-1,1,2,2-tetrachloroethane (98%) was purchased from Alfa-Aesar and was used as received. 1,2-Dibromo-1,1,2,2-tetrafluoroethane (>99.0%, Tokyo Kasei Kogyo Co. Ltd.) is more convenient as a bromination reagent, but is not available in the USA.

8. Addition required approximately 30 sec during which time the temperature rose from -71 °C to -35 °C. If this reagent was added slowly, a significant amount of side products were generated, probably due to self-condensation.

9. The progress of the reaction was monitored by gas chromatography. GC-MS analysis was performed on a JMS-Automass SUN200 instrument equipped with a TC-17 30-m x 0.32-mm x 0.25-µm column under the following conditions: injector temp 250 °C; oven temp 60 °C, 3 min; ramp 8 °C/min; final temp 250 °C; helium gas flow 0.8 mL/min; t_R = 4.56 min (dibromo ester), t_R = 2.13 min (bromo ester). The Checkers did not monitor the reaction, but found that the product could be analyzed by GC-MS (Model HP 5890A GC equipped with an HP 5970 MS detector; column HP-1).

10. A large portion of the aqueous solution froze on contact with the

reaction mixture. The frozen material was allowed to thaw before performing the extraction.

11. A single distillation afforded material that was contaminated by the starting bromo ester (which could be detected by NMR spectroscopy) and the brominating agent (1,2-dibromo-1,1,2,2-tetrachloroethane). The latter impurity was difficult to identify by TLC (not UV active and did not readily stain), GC (did not combust nor did it ionize easily), IR (C-Cl stretches were weak and were obscured by the dibromo ester), or NMR (no protons to observe and the ^{13}C resonance was suppressed by the long relaxation time and the lack of a proton nOe). Combustion analysis proved to be the only reliable method to determine the purity of the distillate. The combustion analysis results clearly demonstrated that the primary distillation did not satisfactorily remove the haloalkane. Both of these impurities were removed by a second distillation (65–67 °C, 9 mmHg) to afford 29.3–29.4 g (74–75%) of analytically pure material. The dibromo ester was a stable compound that can be distilled, purified by silica gel column chromatography, and stored for months in a refrigerator.

12. Physical and spectroscopic properties of **1** are as follows: ^1H -NMR (500 MHz, CDCl_3) δ : 1.36 (t, 3 H, $J = 7.0$), 2.65 (s, 3 H), 4.34 (q, 2 H, $J = 7.0$). ^{13}C -NMR (125 MHz, CDCl_3) δ : 13.7 (q), 37.3, 52.2, 63.7, 166.8. IR (film) cm^{-1} : 1737. EI-MS m/z : 260 ($\text{M}^+ + 2$), 258 (M^+). Anal. Calcd for $\text{C}_5\text{H}_8\text{Br}_2\text{O}_2$: C, 23.10; H, 3.10; Br, 61.48. Found: C, 22.79; H, 2.91; Br, 61.53.

13. A solution of *tert*-butyllithium in pentane was purchased from Sigma Aldrich Company and was titrated by both No-D NMR spectroscopy^{5a} and Gilman double titration.^{5b} This lithium-halogen exchange procedure for the preparation of ynolates⁶ is suitable for a laboratory scale because of its convenience and simplicity. Reductive lithiation with lithium metal catalyzed by naphthalene was better for larger scale.⁷

14. The exothermicity of the reaction subsided substantially after one half of the alkyllithium reagent was added. Thus, the rate of addition was increased after this point without a concomitant increase in internal temperature.

15. Acetyltrimethylsilane (97%) was purchased from Aldrich Chemical

Company and was distilled prior to use (bp 110–116 °C).

16. The reaction of the acylsilane with the lithium ynolate was exothermic and reached a maximum temperature of 39 °C. The addition rate was adjusted to minimize the loss of pentanes from the solution.

17. TLC analysis was performed on Merck silica gel 60 F-254 plates eluting with hexanes/EtOAc, 4:1, R_f = 0.43 (visualized with 254 nm UV lamp).

18. The product was obtained in analytically pure form by distillation on a half scale. The submitters reported that distillation of the product on the full scale led to polymerization.

19. Physical and spectroscopic properties of **2** are as follows: ^1H -NMR (500 MHz, CDCl_3) δ : 0.16 (s, 9 H), 1.90 (q, 3 H, J = 1.0), 1.97 (q, 3 H, J = 1.0), 11.6–12.6 (br s, 1 H). ^{13}C -NMR (125 MHz, CDCl_3) δ : 0.2, 15.4, 20.4, 136.0, 155.6, 175.1. ^{29}Si -NMR (119 MHz, CDCl_3) δ : -5.15. IR (film) cm^{-1} : 2966, 1679. MS (EI) m/z 157 (M^+ – Me). Anal. Calcd. for $\text{C}_8\text{H}_{16}\text{O}_2\text{Si}$ C, 55.77; H 9.36. Found: C, 55.78, H. 9.63.

Safety and Waste Disposal Information


All hazardous materials should be handled and disposed of in accordance with “Prudent Practices in the Laboratory”; National Academy Press; Washington, DC, 1995.

3. Discussion

Ynolates⁸ are ketene anion equivalents and are expected to act as multifunctional carbanions for highly reactive ketenes. Although about 10 reports on the syntheses of ynolates have been published, studies on the versatility of ynolates have remained at a rather basic level, probably due to a lack of convenient and general methods for their synthesis. In *Organic Syntheses*, only one report on ynolates has been published.⁹


The procedure described here demonstrates a practical and convenient method for the generation of ynolates, and a (*Z*)-selective olefination of acylsilanes leading to (*Z*)-2-methyl-3-trimethylsilyl-2-butenic acid.

The method for the generation of ynolates involves thermal cleavage of ester dianions derived from α,α -dibromo esters. We have synthesized aliphatic dibromo esters ($R = \text{Me}, \text{Bu}, i\text{-Pr}, \text{cyclohexyl}, t\text{-Bu}, \text{phenylethyl}, \text{trimethylsilyl}, \text{etc.}$) by the method described herein. Aromatic dibromo esters are synthesized via radical bromination. Two procedures for the generation of ester dianions from the dibromo esters have been developed. Herein, lithium-halogen exchange of dibromo esters via *tert*-butyllithium giving the ester dianions is described. This procedure is very convenient, simple and suitable for laboratory scale. In an alternative procedure, the naphthalene-catalyzed reductive lithiation⁶ of the dibromo esters with lithium metal is employed. As outlined in the following scheme for the generation of ynolates using reductive lithiation catalyzed by naphthalene,¹⁰ the dibromo ester reacts with lithium naphthalenide at $-78\text{ }^{\circ}\text{C}$ to give an α -bromo ester enolate, which is successively lithiated to form the ester dianion. The dianion is cleaved at $0\text{ }^{\circ}\text{C}$ to generate an ynolate by release of lithium ethoxide.




Vinylsilanes are powerful tools in synthetic organic chemistry. Although olefination of acylsilanes is expected to be a useful method for the preparation of vinylsilanes, few reports on this process have been reported due to the unsuitability of conventional olefination reagents. We have developed a method in which olefination of ketones via ynolates provides α,β -unsaturated carboxylates with good to moderate *E/Z* selectivities.¹¹ Using acylsilanes as a ketone in this process, however, extremely high *Z*-selectivity is achieved in the synthesis of various substituted vinylsilanes. The *E/Z* selectivity is controlled in the electrocyclic ring-opening of the intermediate β -lactone enolates, prepared by cycloaddition of carbonyl groups with ynolates. In this step, the transition state would be stabilized by

the orbital interaction between a breaking C-O σ orbital and Si vacant orbitals.¹²



The vinylsilanes produced can be converted into a variety of functionalized multisubstituted alkenes.



1. Institute for Materials Chemistry and Engineering, Kyushu University, 6-1 Kasugakoen, Kasuga, Fukuoka 816-8580, Japan. E-mail: shindo@cm.kyushu-u.ac.jp
2. Institute of Health Biosciences, The University of Tokushima Graduate

School, Sho-machi 1, Tokushima 770-8505, Japan.

3. Shindo, M.; Sato, Y.; Koretsune, R.; Yoshikawa, T.; Matsumoto, K.; Itoh, K.; Shishido, K. *Chem. Pharm. Bull.* **2003**, *51*, 477-478.
4. Shindo, M.; Matsumoto, K.; Mori, S.; Shishido, K. *J. Am. Chem. Soc.* **2002**, *124*, 6840-6841
5. (a) Hoye, T. R.; Aspaas, A. W.; Eklov, B. M.; Voloshin, M. *Org. Lett.* **2004**, *6*, 2567-2570. (b) Gilman, H.; Cartledge, F. K. *J. Organomet. Chem.* **1964**, *2*, 447.
6. Shindo, M.; Sato, Y.; Shishido, K. *Tetrahedron* **1998**, *54*, 2411-2422.
7. Shindo, M.; Koretsune, R.; Yokota, W.; Itoh, K.; Shishido, K. *Tetrahedron Lett.*, **2001**, *42*, 8357-8360.
8. (a) Shindo, M. *Synthesis* **2003**, 2275-2288; (b) Shindo, M. *Tetrahedron* **2007**, *63*, 10-36.
9. Reddy, R.; Kowalski, C. *Org. Syn. Coll. Vol. 9* **1998**, 426-432.
10. Yus, M.; Herrera, R. P.; Guijarro, A. *Chem. Eur. J.* **2002**, *8*, 2574-2584.
11. (a) Shindo, M.; Sato, Y.; Shishido, K. *J. Org. Chem.*, **2000**, *65*, 5443-5445. (b) Shindo, M.; Sato, Y.; Yoshikawa, T.; Koretsune, R.; Shishido, K. *J. Org. Chem.* **2004**, *69*, 3912-3916.
12. For related effects on torquoselectivity, see: (a) Murakami, M.; Miyamoto, Y.; Ito, Y. *Angew. Chem., Int. Ed.* **2001**, *40*, 189-190. (b) Lee, P. S.; Zhang, X.; Houk, K. N. *J. Am. Chem. Soc.* **2003**, *125*, 5072-5079. (c) Mori, S.; Shindo, M. *Org. Lett.* **2004**, *6*, 3945-3948.

Appendix

Chemical Abstracts Nomenclature; (Registry Number)

1,2-Dibromo-1,1,2,2-tetrachloroethane; (630-25-1)

Ethyl 2-bromopropanoate: Propanoic acid, 2-bromo-, ethyl ester: (535-11-5)

Ethyl 2,2-dibromopropionate: Propanoic acid, 2,2-dibromo-, ethyl ester;
(34757-17-0)

Acetyltrimethylsilane; (13411-48-8)



Mitsuru Shindo was born in Tokyo in 1963. He received his PhD degree at University of Tokyo under the supervision of the late Professor K. Koga. From 1990 until 1996, he was an assistant professor at the laboratory of Professor K. Koga. He spent time in Professor R. A. Holton's group (Florida State University) as a postdoctoral fellow from 1992 until 1994. In 1996, he moved to the University of Tokushima as an associate professor working with Professor K. Shishido. In 2005, he joined the faculty at Kyushu University as an associate professor. He is interested in finding new methodologies for synthetic reactions and in the design and synthesis of natural and unnatural bioactive compounds.



Kenji Matsumoto was born in Hiroshima (Japan) in 1978. He graduated from the University of Tokushima in 2001, where he received his Ph.D. in Organic Chemistry in 2006 under the supervision of Professors Mitsuru Shindo and Kozo Shishido. From 2003 until 2006, he was a research fellow of the Japan Society for the Promotion of Science (JSPS). He then worked as a postdoctoral fellow under the direction of Professor Sergey A. Kozmin at the University of Chicago from 2006 to 2007. He is currently an assistant professor at Kyoto University. His research interests are the development of useful synthetic methods and medicinal chemistry.



Kozo Shishido received his Ph.D degree from Tohoku University in 1976 under the direction of the late Professor Tetsuji Kametani. After he was appointed as Assistant Professor at Tohoku University, he had postdoctoral fellowships with Professor A. I. Scott (Texas A&M University) and Professor M. E. Jung (UCLA) from 1978 to 1980. He then moved to University of Tokushima as Associate Professor in 1989. Since 1994, he has been a Full Professor at the Graduate School of Pharmaceutical Sciences, University of Tokushima. His main research interests are devoted to total synthesis of biologically active natural products and development of new preparative methods.





Brian Eklov graduated with a B.S. in Chemistry from the University of Michigan in 1998 and completed his Ph.D. thesis in 2005 under the direction of Professor Thomas. R. Hoyer at the University of Minnesota. His thesis included discussions of a hybrid computational/ ^1H NMR based method for the elucidation of the relative configuration of complex organic molecules, developing the use of non-deuterated solvents for ^1H NMR spectroscopy (No-D NMR), and a total synthesis of (+)-gigantecin, a skipped bis-THF acetogenin. He is presently pursuing post-doctoral research with Professor Scott Denmark at the University of Illinois.

BME 4-075-S-1-C13
Sublimed Material

Pulse Sequence: s2pul

Solvent: CDCl3
Ambient temperature
User: 1-14-87
File: BME_4-075-S-1-C13
INNOVA-500 "sunds1"

Relax. delay 1.000 sec
Pulse 45.0 degrees
Acq. time 1.086 sec
Width 30165.9 Hz
256 repetitions
OBSERVE C13, 125.6487684 MHz
DECOUPLE H1, 499.6981123 MHz
Power 42 dB
continuously on
WALTZ-16 modulated
DATA PROCESSING
Line broadening 1.0 Hz
FT size 65536
Total time 35 min, 47 sec





BME 4-046-d2-2-3
Fractions 2 and 3
from the second distillation

Pulse Sequence: s2pul

Solvent: CDCl₃
Ambient temperature
File: BME_4-046-d2-2-3
INOVA-500 "sunds1"

Pulse 90.0 degrees
Acq. time 4.096 sec
Width 8000.0 Hz
16 repetitions
OBSERVE H1, 500.0729782 MHz
DATA PROCESSING
Line broadening 0.3 Hz
FT size 65536
Total time 1 min, 5 sec



BME 4-046-d2-2-3-C13

Pulse Sequence: s2pul

Solvent: CDCl₃

Ambient temperature

User: 1-14-87

File: BME 4-046-d2-2-3-C13

INOVA-500 "sunds1"

Relax. delay 1.000 sec

Pulse 55.4 degrees

Acq. time 1.024 sec

Width 32000.0 Hz

1024 repetitions

OBSERVE C13, 125.5832523 MHz

DECOUPLE H1, 499.4375375 MHz

Power 44 dB

continuously on

WALTZ-16 modulated

DATA PROCESSING

Line broadening 1.0 Hz

FT size 65536

Total time 34 min, 41 sec

166.769

77.257

77.000



76.743

63.718

52.241

37.256

13.702



BME 4-075-S-1
Sublimed Material

Pulse Sequence: s2pul

Solvent: CDCl₃
Ambient temperature
File: BME 4-075-S-1
INOVA-500 "sunds1"

Pulse 90.0 degrees
Acq. time 4.096 sec
Width 8000.0 Hz
16 repetitions
OBSERVE H1, 500.0729781 MHz
DATA PROCESSING
Line broadening 0.3 Hz
FT size 65536
Total time 1 min, 5 sec

