

TC-6-7-1H

Pulse Sequence: s2pul

Solvent: CDCl<sub>3</sub>

Ambient temperature

File: TC-6-7-1H

INOVA-500 "silicon"

Relax. delay 10.000 sec

Pulse 45.0 degrees

Acq. time 4.096 sec

Width 8000.0 Hz

16 repetitions



OBSERVE H1, 499.4299002 MHz

DATA PROCESSING

Line broadening 0.3 Hz

FT size 65536

Total time 3 min, 45 sec



TC-6-7-13C

Pulse Sequence: s2pul

Solvent: CDCl<sub>3</sub>

Ambient temperature

User: 1-14-87

File: TC-6-7-13C

INOVA-500 "silicon"

Relax. delay 1.000 sec  
Pulse 45.0 degrees  
Acq. time 1.024 sec  
Width 32000.0 Hz  
673 repetitions  
OBSERVE C13, 125.5817522 MHz  
DECOUPLE H1, 499.4315638 MHz  
Power 44 dB  
continuously on  
WALTZ-16 modulated  
DATA PROCESSING  
Line broadening 1.0 Hz  
FT size 65536  
Total time 25 min, 26 sec

