

TC-6-16-1H

Pulse Sequence: s2pul

Solvent: CDCl3

Ambient temperature

File: TC-6-16-1H

INOVA-500 "silicon"

Relax. delay 10.000 sec

Pulse 45.0 degrees

Acq. time 4.096 sec

Width 8000.0 Hz

16 repetitions


OBSERVE H1, 499.4299002 MHz

DATA PROCESSING

Line broadening 0.3 Hz

FT size 65536

Total time 3 min, 45 sec



TC-6-16-1H

Pulse Sequence: s2pul

Solvent: CDCl₃

Ambient temperature

File: TC-6-16-1H

INOVA-500 "silicon"

Relax. delay 10.000 sec

Pulse 45.0 degrees

Acq. time 4.096 sec

Width 8000.0 Hz

16 repetitions

OBSERVE H1, 499.4299002 MHz

DATA PROCESSING

Line broadening 0.3 Hz

FT size 65536

Total time 3 min, 45 sec

